



RGN Antibody, HRP conjugated

Product Code	CSB-PA623001LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15493
Immunogen	Recombinant Human Regucalcin protein (1-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regucalcin (RC) (Gluconolactonase) (GNL) (EC 3.1.1.17) (Senescence marker protein 30) (SMP-30), RGN, SMP30
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.