



AOC3 Antibody, Biotin conjugated

Product Code	CSB-PA624122LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16853
Immunogen	Recombinant Human Membrane primary amine oxidase protein (27-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane primary amine oxidase (EC 1.4.3.21) (Copper amine oxidase) (HPAO) (Semicarbazide-sensitive amine oxidase) (SSAO) (Vascular adhesion protein 1) (VAP-1), AOC3, VAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AOC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.