



PAG 1 Antibody, FITC conjugated

Product Code	CSB-PA644917LC01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29432
Immunogen	Recombinant Bovine Pregnancy-associated glycoprotein 1 protein (54-380AA)
Raised In	Rabbit
Species Reactivity	Bovine
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pregnancy-associated glycoprotein 1 (PAG 1) (EC 3.4.23.-) (Pregnancy-specific protein B) (PSP-B)
Immunogen Species	Bos taurus (Bovine)
Research Area	Others