



SREBF2 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA666241LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q12772 |
| Immunogen | Recombinant Human Sterol regulatory element-binding protein 2 protein (75-164AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sterol regulatory element-binding protein 2 (SREBP-2) (Class D basic helix-loop-helix protein 2) (bHLHd2) (Sterol regulatory element-binding transcription factor 2) [Cleaved into: Processed sterol regulatory element-binding protein 2], SREBF2, BHLHD2 SREBP2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | SREBF2 |