



SREBF2 Antibody, Biotin conjugated

Product Code	CSB-PA666241LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12772
Immunogen	Recombinant Human Sterol regulatory element-binding protein 2 protein (75-164AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sterol regulatory element-binding protein 2 (SREBP-2) (Class D basic helix-loop-helix protein 2) (bHLHd2) (Sterol regulatory element-binding transcription factor 2) [Cleaved into: Processed sterol regulatory element-binding protein 2], SREBF2, BHLHD2 SREBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SREBF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.