



Pollen allergen Phl p 5a Antibody, Biotin conjugated

Product Code	CSB-PA669216HD01EUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q40962
Immunogen	Recombinant Phleum pratense Pollen allergen Phl p 5a protein (1-286AA)
Raised In	Rabbit
Species Reactivity	Phleum pratense
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pollen allergen Phl p 5a, Allergen Phl p Va, Phl p 5a
Immunogen Species	Phleum pratense
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.