



PARP14 Antibody, Biotin conjugated

Product Code	CSB-PA677190LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q460N5
Immunogen	Recombinant Human Poly [ADP-ribose] polymerase 14 protein (1600-1712AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Poly [ADP-ribose] polymerase 14 (PARP-14) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 8) (ARTD8) (B aggressive lymphoma protein 2), PARP14, BAL2 KIAA1268
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PARP14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.