



NEK9 Antibody, HRP conjugated

Product Code	CSB-PA837435LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD19
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek9 protein (821-979AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek9 (EC 2.7.11.1) (Nercc1 kinase) (Never in mitosis A-related kinase 9) (NimA-related protein kinase 9) (NimA-related kinase 8) (Nek8), NEK9, KIAA1995 NEK8 NERCC
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NEK9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.