



NEK9 Antibody, Biotin conjugated

Product Code	CSB-PA837435LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TD19
Immunogen	Recombinant Human Serine/threonine-protein kinase Nek9 protein (821-979AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase Nek9 (EC 2.7.11.1) (Nercc1 kinase) (Never in mitosis A-related kinase 9) (NimA-related protein kinase 9) (NimA-related kinase 8) (Nek8), NEK9, KIAA1995 NEK8 NERCC
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NEK9