



# RAET1E Antibody, HRP conjugated

Product Code	CSB-PA855052LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3ZAR2
Immunogen	Recombinant Human Retinoic acid early transcript 1E protein (31-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAET1E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RAET1E