




SLC30A3 Antibody

Product Code	CSB-PA859940LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99726
Immunogen	Recombinant Human Zinc transporter 3 protein (286-388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 3 (ZnT-3) (Solute carrier family 30 member 3), SLC30A3, ZNT3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC30A3

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA859940LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.