



SLC30A3 Antibody, FITC conjugated

Product Code	CSB-PA859940LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99726
Immunogen	Recombinant Human Zinc transporter 3 protein (286-388AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 3 (ZnT-3) (Solute carrier family 30 member 3), SLC30A3, ZNT3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC30A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.