



SLC30A3 Antibody, Biotin conjugated

Product Code	CSB-PA859940LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99726
Immunogen	Recombinant Human Zinc transporter 3 protein (286-388AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 3 (ZnT-3) (Solute carrier family 30 member 3), SLC30A3, ZNT3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC30A3