



WWOX Antibody, Biotin conjugated

Product Code	CSB-PA873704LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZC7
Immunogen	Recombinant Human WW domain-containing oxidoreductase protein (1-180AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WW domain-containing oxidoreductase (EC 1.1.1.-) (Fragile site FRA16D oxidoreductase) (Short chain dehydrogenase/reductase family 41C member 1), WWOX, FOR SDR41C1 WOX1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WWOX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.