



RPS10P5 Antibody, FITC conjugated

Product Code	CSB-PA878849LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ39
Immunogen	Recombinant Human Putative 40S ribosomal protein S10-like protein (1-176A)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative 40S ribosomal protein S10-like, RPS10P5, RPS10L
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RPS10P5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.