



RBPJL Antibody, Biotin conjugated

Product Code	CSB-PA883361LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBG7
Immunogen	Recombinant Human Recombining binding protein suppressor of hairless-like protein (1-151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Recombining binding protein suppressor of hairless-like protein (Transcription factor RBP-L), RBPJL, RBPL RBPSUHL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBPJL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.