



PRKAB1 Antibody, Biotin conjugated

Product Code	CSB-PA889918LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y478
Immunogen	Recombinant Human 5'-AMP-activated protein kinase subunit beta-1 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5'-AMP-activated protein kinase subunit beta-1 (AMPK subunit beta-1) (AMPKb), PRKAB1, AMPK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKAB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.