



Recombinant Papio hamadryas Apolipoprotein A-I (APOA1)

Product Code	CSB-YP001913ERJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68293
Product Type	Recombinant Protein
Immunogen Species	Papio hamadryas (Hamadryas baboon)
Purity	≥85% (SDS-PAGE)
Sequence	DEPPQT PWDRVKDLVT VYVEALKDSG KDYVSQFEGS ALGKQLNLKL LDNWDSVTST VSKLREQLGP VTQEFWDNLE KETEGLRQEM SKDLEEVKAK VQPYLDDFQK KWQEEMELYR QKVEPLRAEL HEGTRQKLHE LHEKLSPLGE EVRDRARAHV DALRTHLAPY SDELQRQLAA RLEALKENGG ARLAEYHAKA SEHLSTLSEK AKPALEDLRQ GLLPVLESFK VSFLSALEEY TKKLSTQ
Source	Yeast
Target Names	APOA1
Protein Names	Recommended name: Apolipoprotein A-I Short name= Apo-AI Short name= ApoA-I Alternative name(s): Apolipoprotein A1 Cleaved into the following chain: 1. Truncated apolipoprotein A-I
Expression Region	25-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.