




Recombinant Epstein-Barr virus Glycoprotein BDLF3 (BDLF3), partial

Product Code	CSB-EP365881EFA1
Abbreviation	Recombinant Epstein-Barr virus BDLF3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03224
Product Type	Recombinant Protein
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SSTASAGNVTGTTAVTTPSPSASGPSTNQSTTLTTTSAPITTTAILSTNTTTVTF TGTTVTPVPTTSNASTINVTTKVTAQNITATEAGTGTSTGVTSNVVTRSSSTTS ATTRITNATTLAPTLSSKGTSNATKTTAELPTVPDERQPSLSYGLPLWT
Research Area	Others
Source	E.coli
Target Names	BDLF3
Protein Names	Recommended name: Probable membrane antigen GP85
Expression Region	29-186aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	22.7 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

