




mpt64 Antibody

Product Code	CSB-PA14947A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WIN9
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunogenic protein MPT64 (Antigen MPT64), mpt64
Immunogen Species	Mycobacterium tuberculosis
Research Area	Others
Target Names	mpt64

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: mpt64 antibody at 1.5μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 25 kDa

Observed band size: 25 kDa