



Recombinant Human Neuropeptide FF receptor 2 (NPFFR2), partial, Mammalian cell

Product Code	CSB-MP015983HU1
Abbreviation	NPFFR2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0PJM9
Product Type	Recombinant Protein
Purity	≥ 90% as determined by SDS-PAGE.
Research Area	Neuroscience
Source	Mammalian cell
Protein Names	Recommended name: Neuropeptide FF receptor 2 Alternative name(s): G-protein coupled receptor 74 G-protein coupled receptor HLWAR77 Neuropeptide G-protein coupled receptor
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Mol. Weight	50.7kDa
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.