



Recombinant Staphylococcus epidermidis Dihydrofolate reductase (folA)

Product Code	CSB-MP006847SLD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C0P1
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus epidermidis (strain ATCC 12228)
Purity	>85% (SDS-PAGE)
Sequence	MTLSIIVAHD KQRVIGYQNN LPWHLPNDLK HVKQLTTGNT LVMGRKTFNS IGKPLPNRRN VVLTNQASFH HEGVDVINSL DEIKELSGHV FIFGGQTLFE AMIDQVDDMY ITVIDGKFQG DTFPPPYTFE NWEVESSVEG QLDEKNTIPH TFLHLVRRKG K
Source	Mammalian cell
Target Names	folA
Protein Names	Recommended name: Dihydrofolate reductase Short name= DHFR EC= 1.5.1.3
Expression Region	1-161
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.