



# Recombinant Human herpesvirus 1 Uracil-DNA glycosylase (UL2)

<b>Product Code</b>	CSB-MP025641HWY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P10186
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKRACSRSPS PRRRPSSPRR TPPRDGTPPQ KADADDPTPG ASNDASTETR PGSGGEP AAC RSSGPAALLA ALEAGPAGVT FSSSAPPDPP MDLTNGGVSP AATSAPLDWT TFRRVFLIDD AWRPLMEPEL ANPLTAHLLA EYNRRRCQTEE VLPPREDVFS WTRYCTPDEV RVVIIGQDPY HHPGQAHGLA FSVRANVPPP PSLRNVLA AV KNCYPEARMS GHGCLEKWAR DGVLLNNTTL TVKRGAAASH SRIGWDRFVG GVIRRLAARR PGLVFM LWGT HAQN AIRPDP RVHCVLKFSH PSPLSKVPFG TCQHFLVANR YLETRSISPI DWSV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UL2
<b>Protein Names</b>	Recommended name: Uracil-DNA glycosylase Short name= UDG EC= 3.2.2.27 Alternative name(s): UNG
<b>Expression Region</b>	1-334
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.