



# Recombinant Rat Involucrin (Ivl)

<b>Product Code</b>	CSB-BP011922RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P48998
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSHQHTVPVT VPAVVQESLK TVCSPAQTQQ EQTKQPTYP AQCQVFTDTQ EKGFPKHEEK EANPVKDLPE QESEHHQQPG PQKQQLQVKK PEQELQEDEL HSEKQPQEPQ GLLCLGQQQQ REPQEQEQHL RQHQQPQQES QGQGLCLGQQ QDVLAPQELH MGQHQKEKLQ EPELPLGQQQ KTPEEQELIL GEKQQKLHLV ERHQEPQEQE LHHGQKQKQQ QPQEQELQLV QHQQKQHEP ELCLRKQQQQ ESHRELHLG KQQQQESHEP ELHLGKQQHQ ESHEPELHLG KQQHQESCEP ELHLGEQQHQ EQQQHQESCE PELHLGKQQH QETQESELQL GKQQKPHEPD MVLDPKEKQK LHDPELHLGK QHQEQESE LQVEKKQHEK SPEPELHLGK QQELHEPDMT EDQKEKQSLH EPELHLGKQQ ESHEPDMTED QKEKQSLYEP ELHLGKQQEQ QIEYEGYQRS KSLNQLLKQE KASRGQELDD SHLEQEKELL DQRLDQELVN KDEQLERKKH KLENLTQKEK QIKQLVPSTD RVQETQPIQP VKEDSLTTKK QQHSHEVQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	Ivl
<b>Protein Names</b>	Recommended name: Involucrin
<b>Expression Region</b>	1-568
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Involucrin, a component of the keratinocyte crosslinked envelope, is found in the cytoplasm and crosslinked to membrane proteins by transglutaminase. This gene is mapped to 1q21, among calpactin I light chain, trichohyalin, profillaggrin, loricrin, and calcyclin.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.