



Recombinant Human Paxillin (PXN)

Product Code	CSB-EP019112HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49023
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDDL DALLAD LESTTSHISK RPVFLSEETP YSYPTGNHTY QEIAVPPPVP PPPSSEALNG TILDPLDQWQ PSSSRFIHQQ PQSSSPVYGS SAKTSSVSNP QDSVGS PCSR VGEEEHVYSF PNKQKSAEPS PTVMSTSLGS NLSELDRLLL ELNAVQHNPP GFPADANSS PPLPGALSPL YGVPETNSPL GKGAGPLTKE KPKRNGGRGL EDVRPSVESL LDELESSVPS PVPITVNQG EMSSPQRVTS TQQQTRISAS SATRELDLM ASLSDFKIQQ LEQRADGERC WAAGWPRDGG RSSPGGQDEG GFMAQGKTGS SSPPGGPPKP GSQLDSMLGS LQSDLNKLGV ATVAKGVCGA CKKPIAGQVV TAMGKTWHPE HFVCTHCQEE IGSRNFFERD GQPYCEKDYH NLFSPRCYYC NGPILDKVVT ALDRTWHPEH FFCAQCGAFF GPEGFHEKDG KAYCRKDYFD MFAPKCGGCA RAILENYISA LNTLWHPECF VCRECFTPFV NGSFFEHDGQ PYCEVHYHER RGSLSGCGQK PITGRCITAM AKKFHPEHFV CAFCLKQLNK GTFKEQNDKP YCQNCFLKLF C
Source	E.coli
Target Names	PXN
Protein Names	Recommended name: Paxillin
Expression Region	1-591
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.