



# Recombinant Human Paxillin (PXN)

<b>Product Code</b>	CSB-EP019112HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P49023
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDDL DALLAD LESTTSHISK RPVFLSEETP YSYPTGNHTY QEIAVPPPVP PPPSSEALNG TILDPLDQWQ PSSSRFIHQQ PQSSSPVYGS SAKTSSVSNP QDSVGSPPCSR VGEEEHVYSF PNKQKSAEPS PTVMSTSLGS NLSELDRLLL ELNAVQHNPP GFPADEANSS PPLPGALSPL YGVPETNSPL GKGAGPLTKE KPKRNGGRGL EDVRPSVESL LDELESSVPS PVPITVNQG EMSSPQRVTS TQQQTRISAS SATRELDLDEL ASLSDFKIQQ LEQRADGERC WAAGWPRDGG RSSPPGGQDEG GFMAQGKTGS SSPPGGPPKP GSQLDSMLGS LQSDLNKLGV ATVAKGVCGA CKKPIAGQVV TAMGKTWHPE HFVCTHCQEE IGSRNFFERD GQPYCEKDYH NLFSPRCYYC NGPILDKVVT ALDRTWHPEH FFCAQCGAFF GPEGFHEKDG KAYCRKDYFD MFAPKCGGCA RAILENYISA LNTLWHPECF VCRECFTPFV NGSFFEHDGQ PYCEVHYHER RGSLSGSGCQK PITGRCITAM AKKFHPEHFV CAFCLKQLNK GTFKEQNDKP YCQNCFLKLF C
<b>Source</b>	E.coli
<b>Target Names</b>	PXN
<b>Protein Names</b>	Recommended name: Paxillin
<b>Expression Region</b>	1-591
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.