



Recombinant Lemur catta Involucrin (IVL)

Product Code	CSB-MP011922LOK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14590
Product Type	Recombinant Protein
Immunogen Species	Lemur catta (Ring-tailed lemur)
Purity	>85% (SDS-PAGE)
Sequence	MSQQHTLPVT LPPTLSQELL KNVSPADIQ QEQRKQPTPL PAPCQKVLSE LPVAVPSKHE EKHATPVKGL LEQECGQLQQ QEPQEQEVHL AKHQELQELQ EQELHLGKQP EPQEQELHLG KQQQQQESQE QELYLGKQPE PQDQELHLGK QQAPQEQELH LGKQPEPQEQ ELHLGKQPEP QDQELYLGKR LEPQEQELHL GKQQQQQESQ EQELDLGKQP EPQDQELHLG KQQAPQEQEL HLGKQPEPQE QELHLVKQQE PQDQELHLGK RLEPQEQELH LGKQQQPQEQ KLHPGAAAAA GVTGAGPAAS KAARRATGAG TAPGAAAAA GATGAGTAAT APATAEERQK AESLEQQLEQ EKAQREEQLK EQLEEKKRIL DQQLDQEVAK RYEQLGVKKE QLLQPLGQQE GQLEKPVFVP APGQVQDIQP PQPPKGEVLL PAEQQQEPEV
Source	Mammalian cell
Target Names	IVL
Protein Names	Recommended name: Involucrin
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.