



Recombinant Human Protein ERGIC-53 (LMAN1)

Product Code	CSB-MP012991HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49257
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DGVGGDPAVALPHRRFEYKYSFKGPHLVQSDGTVPFWAHAGNAIPSSDQIRV APSLKSQRGSVWTKTKAAFENWEVEVTFRTGRGRIGADGLAIWYAENQGLE GPVFGSADLWNGVGIFDSDNDGKKNPAIVIIIGNNGQIHVDHQNDGASQAL ASCQRDFRNKPYVRAKITYYQNTLTVMINNGFTPDKNDEYFCAKVENMIIPAQ GHFGISAATGGLADDHDVLSFLTFQLTEPGKEPPTPDKEISEKEKEKYQEEFEH FQQELDKKKEEFQKGHPDLQGQPAEEIFESVGDRELRLQVFEGQNRHLEIKQL NRQLDMILDEQRRYVSSLTEEISKRGAGMPGQHGQITQQELDTVVKTQHEILR QVNEMKNSMSETVRLVSGMQHPGSAGGVYETTQHFIDIKEHLHIVKRDIDNLV QRNMPSNEKPKCELPFPSCLS
Source	Mammalian cell
Target Names	LMAN1
Protein Names	Recommended name: Protein ERGIC-53 Alternative name(s): ER-Golgi intermediate compartment 53 kDa protein Gp58 Intracellular mannose-specific lectin MR60 Lectin mannose-binding 1
Expression Region	31-477aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Target Details	This protein is a type I integral membrane protein localized in the intermediate region between the endoplasmic reticulum and the Golgi, presumably recycling between the two compartments. The protein is a mannose-specific lectin and is a member of a novel family of plant lectin homologs in the secretory pathway of animal cells. Mutations in the gene are associated with a coagulation defect. Using positional cloning, the gene was identified as the disease gene leading to combined factor V-factor VIII deficiency, a rare, autosomal recessive disorder in which both coagulation factors V and VIII are diminished.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.