



Recombinant *Cricetulus griseus* Glyceraldehyde-3-phosphate dehydrogenase (GAPDH), partial

Product Code	CSB-MP009232DXU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17244
Product Type	Recombinant Protein
Immunogen Species	<i>Cricetulus griseus</i> (Chinese hamster) (<i>Cricetulus barabensis griseus</i>)
Purity	>85% (SDS-PAGE)
Sequence	VKVGVNGFG RIGRLVTRAA FTSGKVEVVA INDPFIDLNY MVYMFQYDST HGKFKGTVKA ENGLVINGK AITIFQERDP ANIKWGDAGA EYVVESTGVF TTMEKAG AHL KGGAKRVIIS APSADAPMFV MGVNQDKYDN SLKIVSNASC TTNCLAPLAK VIHDNFGIVE GLMTTVHAIT ATQKTVDGPS GKLWRDGRGA AQNIIPASTG AAKAVGKVIP ELNGKLTGMA FRVPTPNVSV VDLTCRLEKP AKYEDIKKVV KQASEGPLKG ILGYTEDQVV SCDFNSDSHS STFDAGAGIA LNDNFVKLIS WYDNEFGYSN RVVDLMAYMA SKE
Source	Mammalian cell
Target Names	GAPDH
Protein Names	Recommended name: Glyceraldehyde-3-phosphate dehydrogenase Short name= GAPDH EC= 1.2.1.12Alternative name(s): Peptidyl-cysteine S-nitrosylase GAPDH EC= 2.6.99.-
Expression Region	2-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.