



# Recombinant *Macropus eugenii* Phosphoglycerate kinase 1 (PGK1)

<b>Product Code</b>	CSB-MP017856MPE
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P29408
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macropus eugenii</i> (Tamar wallaby) ( <i>Notamacropus eugenii</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SLSKKLTL D KVDLKGK RVI MRVDFNV PMK NKEITNNQRI KAALPSIN YC LDNGAKSVVL MSHLGRPDGV PMPDKYSLEP VAAELK SLLG KDVLFLKDCV GPEVEQACAN PPNGSIILLE NLRFHVEEEG KGKDASGNKI KAEPAKMEAF QASLSKLG DV YVNDAFGTAH RAHSSMVG VN LPQKACGFLM KKELTYFAKA LDSPERPFLA ILGGAKVADK IQLINMLDK VNEMIIGGGM GFTFLKVLNN MEIGTSLFDE EGAKIVK DLM AKA EKNGVKI TLPVDFVTAD KFDENAKTGQ ATLASGIPAG WMGLDCGPKS SKKYVEVVTW AKQIVWNGPV GVFEWEEFAR GTKELMNNVV EATKRCITI IGGGDTATCC AKWNTEDKVS HVSTGGGASL ELLEGKVLPG VSTLNNV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PGK1
<b>Protein Names</b>	Recommended name: Phosphoglycerate kinase 1 EC= 2.7.2.3
<b>Expression Region</b>	2-417
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.