



Recombinant *Xenopus laevis* Follistatin (fst)

Product Code	CSB-MP009024XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31515
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	G NCWLQQSKNG RCQVLYRTEL SKEECCKTGR LGTSWTEEDV PNSTLFKMWI FHGGAPHCIP CKETCENVDC GPGKKCKMKNK KNKPRCVCAP DCSNITWKGS VCGIDGKTYK DECALLKAKC KGVPELDVQY QGKCKKTCRD VLCPGSSSCV VDQTNNAYCV TCNRICPEPT SPDQYLCGND GITYGSACHL RKATCLLGRS IGLAYEGKCI KAKSCEDIQC SAGKKCLWDS RVGRGRCALC DDLCGESKSD DTVCASDNTT YPSECAMKQA ACSTGILLEV KHSGSCNSIV EDTEEEEEEE EPDYSFVISS W
Source	Mammalian cell
Target Names	fst
Protein Names	Recommended name: Follistatin Short name= FS Alternative name(s): Activin-binding protein XFS-319
Expression Region	30-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.