



Recombinant *Saccharomyces cerevisiae* Uroporphyrinogen decarboxylase (HEM12)

Product Code	CSB-EP025678SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P32347
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGNFPAPKND LILRAAKGEK VERPPCWIMR QAGRYLPEYH EVKNNRDFEQ TCRDAEIASI ITIQPVRRYR GLIDAAIIFS DILVIPQAMG MRVEMLEGKG PHFPEPLRNP EDLQTVLDYK VDVLKELDWA FKAITMTRIK LDGEVPLFGF CGGPWTLMVY MTEGGGSRFL RFAKQWINMY PELSHKLLQK ITDVAVEFLS QQVVAGAQL QVFESWGGEL SSVDFDEFSL PYLRQIAERV PKRLQELGIM EQIPMIVFAK GSWYALDKLC CSGFDVVSLD WSWDPREAVK INKNRVTLQG NLDPGVMYGS KEVITKKVKQ MIEAFGGGKS RYIVNFGHGT HPFMDPDVIK FFLEECHRIG SK
Source	E.coli
Target Names	HEM12
Protein Names	Recommended name: Uroporphyrinogen decarboxylase Short name= UPD Short name= URO-D EC= 4.1.1.37
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.