



Recombinant *Saccharomyces cerevisiae* Glutathione peroxidase-like peroxiredoxin 1 (GPX1)

Product Code	CSB-BP009866SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36014
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MQEFYSFSPI DENGNPFFFN SLRNKVV LIV NVASHCAFTP QYKELEYLYE KYKSHGLVIV AFPCGQFGNQ EFEKDKEINK FCQDKYGVTF PILHKIRCNG QKQDPVYKFL KNSVSGKSGI KMIKWNFEKF VVDRNGKVVK RFSCMTRPLE LCPIIEELLN QPPEEQI
Source	Baculovirus
Target Names	GPX1
Protein Names	Recommended name: Glutathione peroxidase 1 EC= 1.11.1.9
Expression Region	1-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.