



Recombinant Human 1-acyl-sn-glycerol-3-phosphate acyltransferase beta (AGPAT2), partial, Baculovirus

Product Code	CSB-BP001450HU1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O15120
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	AGPAT2
Protein Names	Recommended name: 1-acyl-sn-glycerol-3-phosphate acyltransferase beta EC= 2.3.1.51 Alternative name(s): 1-acylglycerol-3-phosphate O-acyltransferase 2 Short name= 1-AGP acyltransferase 2 Short name= 1-AGPAT 2 Lysophosphatidic aci
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.