




Recombinant Mouse Dihydroorotate dehydrogenase (quinone), mitochondrial (Dhodh), partial

Product Code	CSB-EP006852MO1
Relevance	Catalyzes the conversion of dihydroorotate to orotate with quinone as electron acceptor.
Abbreviation	Recombinant Mouse Dhodh protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O35435
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	TGDDHFYAEYLMPALQRLLDPESAHLRAVRVISLGLLPRATFQDSNMLEVRVL GHKFRNPVGIAAGFDKHGEAVDGLYKLGFGFVEVGSVTPQPQEGNPRPRVFR LPEDQAVINRYGFNSHGLSAVEHRLRARQQKQTQLTTDGLPLGINLGKKNKTSV DAAADYVEGVRILGPLADYLVVNVSSPNTAGLRSLQGKTELRRLLSKVLQERD ALKGPQKPAVLVKIAPDLTAQDKEDIASVARELGIDGLIITNTTVSRPVGLQGAL RSETGGLSGKPLRDLSTQTIREMYALTQGTIPIIGVGGVSSGQDALEKIQAGAS LVQLYTALTFLGPPVVARVKRELEALLKERGFNTVTD AIGVDHRR
Research Area	Metabolism
Source	E.coli
Target Names	Dhodh
Protein Names	Dihydroorotate oxidase
Expression Region	31-395aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	45.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.