



Recombinant Schizosaccharomyces pombe

Argininosuccinate lyase (arg7)

Product Code	CSB-BP002213SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40369
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAEKSSKKLW GGRFSGATDP LMAEFNKSIY SGKEMCEEDV IGSMAYAKAL CQKNVISEEE LNSILKGLEQ IQREWNSGQF VLEPSDEDVH TANERRLTEI IGDVAGKLHT GRSRNDQVTT DLRLWLCKRI KEVEVYVINL LKVFTNRAEM EIDVIMSGYT HLQRAQPVRW SHFLMSHALP LLGDLGRLRQ LYTRVSQLPL GAGALAGNPF NVDREFLRKE LGFEGIIMNS MNAVGD RDFV IEFMFWAGMV MLHISRFAED LIYSSSEFG FVTLS DAYST GSSIMPQKKN PDSLELLRGK SGRVLGDMIG LMITVKG TPT TYNKDLQEDK EPLFDAFKTV SDSLQILTGV VSTLTINPTK IAESLTPDLL ATDLAEYLVR KGLPFRQTHH ISGSAVRMAE ERNTTLDKLS VSDLQSLHPL FDEDVSKVFN YEESVEKRCS IGGTAKHCVQ EQIEHIRSAI L
Source	Baculovirus
Target Names	arg7
Protein Names	Recommended name: Argininosuccinate lyase Short name= ASAL EC= 4.3.2.1 Alternative name(s): Arginosuccinase
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.