



Recombinant Octodon degus Glucagon (GCG)

Product Code	CSB-EP009315ODV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P22890
Product Type	Recombinant Protein
Immunogen Species	Octodon degus (Degu) (Sciurus degus)
Purity	>85% (SDS-PAGE)
Sequence	HPLQDTEEKP RSFSTSQTDL LDDPDQMNEK KRHSQGTFTS DYSKFLDTRR AQDFLDWLKN TKRNRNEIA
Source	E.coli
Target Names	GCG
Protein Names	Recommended name: Glucagon Cleaved into the following 8 chains: 1. Glicentin 2. Glicentin-related polypeptide Short name= 3. GRPP 4. Oxyntomodulin Short name= 5. OXM Short name= 6. OXY 7. Glucagon 8. Glucagon-I
Expression Region	21-89
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.