



Recombinant *Drosophila melanogaster* Calreticulin (Crc)

Product Code	CSB-YP004458DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29413
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	SAE VYLKENFDNE NWEDTWIYSK HPGKEFGK FV LTPGTFYNDA EADKGIQTSQ DARFYAASRK FDGFSNEDKP LVVQFSVKHE QNIDCGGGYV KLFDCSLDQT DMHGESPYEI MFGPDCGPG TTKVHVIFSY KGKNHLISKD IRCKDDVYTH FYTLIVRPDN TYEVLIDNEK VESGNLEDDW DFLAPKKIKD PTATKPEDWD DRATIPDPDD KKPEDWDKPE HIPDPDATKP EDWDDemdGE WEPPMIDNPE FKGEWQPKQL DNPNYKGAW E HPEIANPEYV PDDKLYLRKE ICTLGFDLWQ VKSGTIFDNV LITDDVELAA KAAAEVKNTQ AGEKMKKEAQ DEVQRKKDEE EAKKASDKDD EDEDDDDDEEK DDESKQDKDQ SEHDEL
Source	Yeast
Target Names	Crc
Protein Names	Recommended name: Calreticulin
Expression Region	18-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.