



Recombinant Anser anser anser Fatty acid synthase (FASN)

Product Code	CSB-BP008435AKX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36189
Product Type	Recombinant Protein
Immunogen Species	Anser anser anser (Western greylag goose)
Purity	>85% (SDS-PAGE)
Sequence	MEDVVIAGIA GKLPESENLQ EFWEKLLNGV DMVTEDDRRW KPGMYGLPKR NGKLDISKF DASFFGVHPK QAHTMDPQLR LLEVSYEAI LDGGINPATL RGTDTGVWVG ASGSEAGEAF SQDPEQLLGY SMIGCQRAMF ANRISYFYDF KGPSSLIDTA CSSSLMALEN AYKAIRNGRC SAAVVGGVNL LKPNNTSVQF MKLGMLSPDG ACKVFDASGD GYCRSEAVVV VLLTKKSMK RIYATIVNAG SNTDGFKEQG VTFPSGDMQR QLVSSLHREC GIKPGDIEYV ETHGTGTVKVG DPQEVNGLAD LFCQCEREPL LIGSTKSNMG HPEPASGLAA LAKVVLSLEH GL
Source	Baculovirus
Target Names	FASN
Protein Names	Recommended name: Fatty acid synthase EC= 2.3.1.85
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.