



Recombinant Mouse Homeobox even-skipped homolog protein 2 (Evx2)

| | |
|--------------------------|--|
| Product Code | CSB-EP007874MO-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P49749 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MMERIRKEMI LMERGLHSPT AGKRFSSLSD SAGGAVLEAL ENSQHPARLS PRLPSAPLHG ALGDLPAK GK FEIDTLFNLQ HPSSESTVSS EIASATESRK KPSHYSEAAA EADMSSDVEV GCSALRSPSG LGAAPLKENN AKGYTESGSV AGTTTSASGS GLGSLHGGGG GGN SGAAALG GSGSGSGADQ VRRYRTAFTR EQIARLEKEF YRENYVSRPR RCELAALNL PETTIKWWFQ NRRMKDKRQR LAMSWPHPAD PSFYTYMMTH AAATGSLPYP FHSVPLHYY PHVGVTA AAAASGAAAA ASSPFATSIR PLDTFRALSH PYSRPELLCS FRHPGLYQAP AAAAGLNSAA SAAAAAAAAA AAASSAAAAG APPSGSSAPC SCLSCHSSQS AAAAAAAAAA ALGSRGGGGS GGGGGGGAGT AGGSDFGCSA AAPRSESGFL PYSAAVLSKT AVSPDQRDE APLTR |
| Source | E.coli |
| Target Names | Evx2 |
| Protein Names | Recommended name: Homeobox even-skipped homolog protein 2 Alternative name(s): EVX-2 |
| Expression Region | 1-475 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |