



Recombinant Mouse Homeobox even-skipped homolog protein 2 (Evx2)

Product Code	CSB-MP007874MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49749
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MMERIRKEMI LMERGLHSPT AGKRFSSLSD SAGGAVLEAL ENSQHPARLS PRLPSAPLHG ALGDLPAKGK FEIDTLFNLQ HPSSESTVSS EIASATESRK KPSHYSEAAA EADMSSDVEV GCSALRSPSG LGAAPLKENN AKGYTESGSV AGTTTSASGS GLGSLHGGGG GGN SGAAALG GSGSGSGADQ VRRYRTAFTR EQIARLEKEF YRENYVSRPR RCELAALNL PETTIKWWFQ NRRMKDKRQR LAMSWPHPAD PSFYTYMMTH AAATGSLPYP FHSVPLHYY PHVGVTA AAAASGAAAA ASSPFATSIR PLDTFRALSH PYSRPELLCS FRHPGLYQAP AAAAGLNSAA SAAAAAAAAA AAASSAAAAG APPSGSSAPC SCLSCHSSQS AAAAAAAAAA ALGSRGGGGS GGGGGGGAGT AGGSDFGCSA AAPRSESGFL PYSAAVLSKT AVSPDQRDE APLTR
Source	Mammalian cell
Target Names	Evx2
Protein Names	Recommended name: Homeobox even-skipped homolog protein 2 Alternative name(s): EVX-2
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.