



# Recombinant Human 5,6-dihydroxyindole-2-carboxylic acid oxidase (TYRP1), partial, Baculovirus

<b>Product Code</b>	CSB-BP025399HU1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P17643
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Source</b>	Baculovirus
<b>Target Names</b>	TYRP1
<b>Protein Names</b>	Recommended name: 5,6-dihydroxyindole-2-carboxylic acid oxidase Short name= DHICA oxidase EC= 1.14.18.-Alternative name(s): Catalase B Glycoprotein 75 Melanoma antigen gp75 Tyrosinase-related protein 1 Short name= TRP S
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
<b>Protein Length</b>	Partial
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.