



Recombinant Mouse B2 bradykinin receptor (Bdkrb2), partial, Escherichia Coli

| | |
|--------------------------|---|
| Product Code | CSB-EP002652MO1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P32299 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Source | E.coli |
| Target Names | Bdkrb2 |
| Protein Names | Recommended name: B2 bradykinin receptor Short name= B2R Short name= BK-2 receptor |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-SUMO-tagged and C-terminal Myc-tagged N-terminal His-B2M-tagged and C-terminal Myc-tagged N-terminal His-GST-tagged and C-terminal Myc-tagged N-terminal His-tagged N-terminal His-SUMO-tagged N-terminal His-B2M-tagged N-terminal His-GST-tagged N-terminal GST-tagged Tag-Free |
| Protein Length | Partial |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |