



Recombinant Human Gamma-glutamyltranspeptidase 2 (GGT2), partial, Mammalian cell

Product Code	CSB-MP009395HU1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36268
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	GGT2
Protein Names	Recommended name: Gamma-glutamyltranspeptidase 2 Short name= GGT 2 EC= 2.3.2.2 Alternative name(s): Gamma-glutamyltransferase 2 Glutathione hydrolase 2 EC= 3.4.19.13 Cleaved into the following 2 chains: 1. Gamma-glutamyltranspe
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.