



Recombinant Human Neuropeptide Y receptor type 1 (NPY1R), partial, Mammalian cell

| | |
|--------------------------|---|
| Product Code | CSB-MP016035HU1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P25929 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Source | Mammalian cell |
| Target Names | NPY1R |
| Protein Names | Recommended name: Neuropeptide Y receptor type 1 Short name= NPY1-R |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free |
| Protein Length | Partial |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |