



Recombinant Human Dipeptidyl peptidase 4 (DPP4), partial, Escherichia Coli

Product Code	CSB-EP007139HU1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27487
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Source	E.coli
Target Names	DPP4
Protein Names	Recommended name: Dipeptidyl peptidase 4 EC= 3.4.14.5 Alternative name(s): ADABP Adenosine deaminase complexing protein 2 Short name= ADCP-2 Dipeptidyl peptidase IV Short name= DPP IV T-cell activation antigen CD26 TP103
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-SUMO-tagged and C-terminal Myc-tagged N-terminal His-B2M-tagged and C-terminal Myc-tagged N-terminal His-GST-tagged and C-terminal Myc-tagged N-terminal His-tagged N-terminal His-SUMO-tagged N-terminal His-B2M-tagged N-terminal His-GST-tagged N-terminal GST-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.