



Recombinant Rat Alpha-2A adrenergic receptor (Adra2a), partial, Mammalian cell

Product Code	CSB-MP001388RA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P22909
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	Adra2a
Protein Names	Recommended name: Alpha-2A adrenergic receptor Alternative name(s): Alpha-2A adrenoreceptor Short name= Alpha-2A adrenoceptor Short name= Alpha-2AAR Alpha-2D adrenergic receptor CA2-47
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.