



Recombinant Mouse Matrix metalloproteinase-14 (Mmp14), partial, Baculovirus

Product Code	CSB-BP014661MO1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53690
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	Mmp14
Protein Names	Recommended name: Matrix metalloproteinase-14 Short name= MMP-14 EC= 3.4.24.80 Alternative name(s): MMP-X1 MT-MMP Membrane-type matrix metalloproteinase 1 Short name= MT-MMP 1 Short name= MTMMP1 Membrane-type-1 matrix
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.