



Recombinant Chicken Frizzled-7 (FZD7), partial, Mammalian cell

| | |
|--------------------------|---|
| Product Code | CSB-MP009110CH1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | O57329 |
| Product Type | Recombinant Protein |
| Immunogen Species | Gallus gallus (Chicken) |
| Purity | >85% (SDS-PAGE) |
| Source | Mammalian cell |
| Target Names | FZD7 |
| Protein Names | Recommended name: Frizzled-7 Short name= Fz-7 Short name= cFz-7 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free |
| Protein Length | Partial |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |