



Recombinant Enteroctopus dofleini Rhodopsin (RHO), partial, Baculovirus

Product Code	CSB-BP019681ELZ1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P09241
Product Type	Recombinant Protein
Immunogen Species	Enteroctopus dofleini (North Pacific giant octopus) (Octopus dofleini)
Purity	>85% (SDS-PAGE)
Source	Baculovirus
Target Names	RHO
Protein Names	Recommended name: Rhodopsin
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.